

# nano3 *Pediatrics*





## nano3 *Pediatrics*

We have introduced the world of sound to the individuals that have profound hearing loss and make their lives easier with the slogan of '**more powerful than the most powerful**'.

We continue to serve worldwide with the most reliable and distinctive software.

### For Those Who Want Special Devices with Unique Technology

**We produce the technologies that give meaning to sounds**

nano3 *Pediatrics* we developed for babies and children offers variety of sounds to the little users and introduces existence of the sound and its impacts to them. It contributes to the development of mind and formation of perception and speaking abilities.

## Special For Our Little Users

**nano3** *Pediatrics* we developed has been designed in the way that the children can get maximum benefit from the technology and get solutions for their needs to build a perfect future.

## Some Special Basic Features

### Special And Distinctive Feedback Manager

The ear canal of babies and children are narrower, shorter and straighter than the ear canals of adults. To make mold for this type of ear canal is very hard and does not give good results for most of the time. The molds that do not fit in the ear canal properly causes feedback especially for powerful hearing aids.

We developed algorithms that reduce the impact of the problem and even remove the problem itself completely most of the time.

### Treble Sounds Can Be Heard

Instead of seeking additional and complicated technologies to hear treble sounds, it is possible to hear the treble sounds by **nano3** technology very well without reducing gain.

Listening to music and hearing the sound of birds are no longer a dream.

### Maximum Gain

Provides a balance between frequencies and strengthens required speech frequencies.

## Powerful Performance

Increases the pleasure of hearing adapting to the environment and makes it comfortable to use the hearing aid.

## More Clear Voices


Our new generation technology can detect and balance the needed signals. It recognizes and amplifies frequency based weak toned syllables.

## Difficult Acoustic Conditions

Creates a comfortable environment with the sound suppressor management in difficult acoustic places. Suppresses the noise and amplifies the voice to the level where it can be understood.

## Stable and Durable

nano3 *Pediatrics* hearing aids have been developed for children. It is more resistant to water, dust, humidity and the situations where the hearing aid falls down.





## Battery Consumption

nano3 *Pediatrics* provides the joy of long lasting hearing with low battery consumption. The user keeps enjoying hearing with low cost without changing the battery for a long time.

## Battery Light Signal

A button which lights every five seconds is placed on the hearing aid assuming that children cannot realize whether the battery is run out unlike the other hearing aids. It helps parents check on the battery consumption.

Besides, a speaking technology that gives a warning about battery's running out is integrated.



# Distinctive Models for Pediatrics

There are different types of devices depending on the hearing loss. Hearing specialists suggest the appropriate device.

## Size And Type

New models are more unique, more elegant and more powerful. There are different options for moderate, severe and profound hearing loss.



### SPn

For moderate and severe hearing loss



### HPn

For severe and profound hearing loss



### HPn-L

For severe and profound hearing loss

Head Office | Manufactory

EU



Caliislar Cad. Hanimeli Sokak, No.1/10  
34180 Bahcelievler-Istanbul/Turkey  
Tel.+90 212 557 91 68



Karamehmet Mh. Adnan Arisoy Bulv.  
European Free Trade Zone, No.13/305  
Ergene-Tekirdag/Turkey  
[info@eartech.com](mailto:info@eartech.com) | [eartech.com](http://eartech.com)



Mühlengrund 2, 45355  
Essen / Germany  
[kontakt@eartech.com](mailto:kontakt@eartech.com)  
[eartech.com](http://eartech.com)